

AmonMed Covid-19 Antigen Lollitest

Laientest / Selbsttest mit CE1434 (1er verpackt)

Omikron-Variante (B.1.1.529) erkennen

Hersteller	Xiamen AmonMed Biotechnology Co., Ltd
Rep	SUNGO Europe B.V.
CE	CE1434 seit 15.10.2021
BfArM Nummer	AT1279/21
Paul-Ehrlich-Institut	evaluiert
EU List	Device #1763
HSC Common List	ja
Sensitivität Cq \leq 25	100%
Sensitivität Cq 25-30	87%
Sensitivität \geq 30	30%
Gesamtsensitivität	80%

Varianten (SKU)	1er verpackt
Inhalt pro Karton/VPE	500 St.
Abmessungen Karton	
Gebrauchsanleitung	auf Deutsch



AmonMed™

COVID-19 Antigen Rapid Test Kit (Colloidal Gold)

Immunochromatographic selftest for the rapid qualitative detection of SARS-CoV-2 antigen in human saliva swab.

DE: Ein immunochromatographischer Selbsttest für den schnellen qualitativen Nachweis von SARS-CoV-2-Antigen im menschlichen Speichel.

PL: Immunochromatograficzny test do samokontroli w celu szybkiego, jakościowego wykrywania antygenu SARS-CoV-2 w wymazie z ludzkiej śliny

Distributed by
 **CE** 1434 **IVD**   21°C-30°C  1 

Xiamen AmonMed Biotechnology Co., Ltd.
 Unit 503, 120 Xinyuan Road, Haicang District, Xiamen, Fujian, China
 Tel:+86-592-6519307 www.amonmed.com

EC REP **SUNGO Europe B.V**
 Olympisch Stadion 24, 1076DE
 Amsterdam, Netherlands



Antigen-Tests auf SARS-CoV-2 zur Eigenanwendung

die Gegenstand des Anspruchs nach §1 Satz 1 Coronavirus-Testverordnung (TestV) sind („Selbsttests“)

Suchen: Alle Textspalten		Los	Aktionen	Hersteller		Europäischer Bevollmächtigter		
Test-ID	Name des Tests	Evaluierung PEI	Omikron-Erkennung entsprechend der Bridging-Prüfung des PEI	Name ↓	Land	Name	Land	Probennahme
AT1279/21	COVID-19 Antigen Rapid Test Kit (Colloidal Gold)	Ja	Ja	Xiamen AmonMed Biotechnology Co., Ltd.	CN	SUNGO Europe B.V.	NL	Speichel

AT-Nr. / AT-No.	AT-Nr. Selbsttest / AT-No. self-test	Ref-Nr. / ID-No. *	Hersteller / Manufacturer	Testname / Test name	Zielantigen / target antigen	Cq ≤25	Cq 25-30	Cq ≥30	Gesamt-Sensitivität / total sensitivity
AT211/20		n.a.	Shenzhen Zhenrui Biotech Co., Ltd.	Zhenrui COVID-19 (SARS-COV-2) Antigen Test Kits	N	82,4%	13,0%	0,0%	34,0%
AT941/21		T3701W, T3702W, T3706W	Sichuan Xincheng Biological Co., Ltd.	SARS-CoV-2 Antigen Assay Kit by Latex Immunochromatography method	N	100,0%	95,0%	40,0%	86,0%
AT010/20		CAGT025E0	Sugentech, Inc.	SGT-flex COVID-19 Ag	N	100,0%	73,9%	0,0%	68,0%
AT789/21		COVID19AGVCG	SureScreen Diagnostics Ltd.	COVID-19 Antigen Rapid Test Cassette (Nasal Swab)	N	100,0%	20,0%	0,0%	48,0%
AT765/21		n.a.	Surge Medical Inc.	COVID-19 Antigen Test Kit	N	100,0%	40,0%	0,0%	56,0%
AT830/21		n.a.	Suzhou Soochow University Saier Immuno Biotech Co., Ltd.	InstantSure Covid-19 Ag CARD	N	100,0%	100,0%	80,0%	96,0%
AT848/21		04A024	TBG Biotechnology Xiamen Inc.	TBG SARS-CoV-2 Antigen Rapid Test	N	95,0%	80,0%	10,0%	72,0%
AT379/20		2276-20	Toda Pharma	Toda Coronadiag Ag	N	100,0%	95,7%	40,0%	86,0%
AT044/21		C011906	Triplex International Biosciences (China) Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Kit	N	100,0%	87,0%	20,0%	78,0%
AT492/20	AT1222/21	RTC03	TürkLab Tibbi Malzemeleri San. Ve Tic. A.S	Rapidan Tester COVID-19 Ag Test	N	90,0%	10,0%	0,0%	40,0%
AT075/21		012G521	utili med Products (Deutschland) GmbH	COVID-19 Antigen Speichelttest (Immunochromatographie)	N	100,0%	95,7%	20,0%	82,0%
AT507/20		VSCD02	Vitrosens Biyoteknoloji Ltd. Sti	RapidFor SARS-CoV-2 Rapid Antigen Test Colloidal Gold	N	100,0%	30,4%	0,0%	48,0%
AT711/21		n.a.	Weihai Kangzhou Biotechnology Engineering Co., Ltd.	Kanzone COVID-19 Antigen Rapid Test	N	90,0%	5,0%	0,0%	38,0%
AT430/20		W-Ag03-20	Wuhan EasyDiagnosis Biomedicine Co., Ltd.	COVID-19 (SARS-CoV-2) Antigen Test Kit	N	100,0%	73,9%	0,0%	68,0%
AT997/21		FP-318	Wuhan HealthCare Biotechnology Co., Ltd.	SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	N	90,0%	0,0%	0,0%	36,0%
AT172/20		SF24025	Wuhan Life Origin Biotech Joint Stock Co., Ltd.	SARS-CoV-2 Antigen Assay Kit (Immunochromatography)	n.a.	100,0%	56,5%	0,0%	60,0%
AT429/20		CoV2Ag-25	Wuhan UNscience Biotechnology Co., Ltd.	SARS-CoV-2 Antigen Rapid Test Kit	N	88,2%	17,4%	0,0%	38,0%
AT847/21		COV-S31	Wuxi Biohermes Bio & Medical Technology Co., Ltd.	SARS-CoV-2 Antigen Test Kit (Lateral Flow Assay)	N	100,0%	60,0%	0,0%	64,0%
AT246/21	AT1279/21	CG01Ag-25	Xiamen AmonMed Biotechnology Co., Ltd.	COVID-19 Antigen Rapid Test Kit (Colloidal Gold)	N	100,0%	87,0%	30,0%	80,0%



CERTIFICATE

EC Certificate No. 1434-IVDD-467/2021

**EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices**

Polish Centre for Testing and Certification certifies
that manufactured by:

Xiamen AmonMed Biotechnology Co., Ltd
Unit 503, 120 Xinyuan Road, Haicang District,
Xiamen, Fujian, China

in vitro diagnostic medical devices
for self-testing

COVID-19 Antigen Rapid Test Kit (Colloidal Gold)
Saliva specimen

CG01Ag-01S-ST, CG01Ag-05S-ST, CG01Ag-25S-ST

in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 15.10.2021 to 27.05.2024

The date of issue of the Certificate: 15.10.2021

The date of the first issue of the Certificate: 15.10.2021



Issued under the Contract No. MD-128/2021
Application No: 233/2021
Certificate bears the qualified signature.
Warsaw, 15/10/2021
Module A1
FBM-30-E_10

Anna
Małgorzata
Wyroba
Elektronicznie
podpisany przez Anna
Małgorzata Wyroba
Data: 2021.10.15
13:27:16 +02'00'
Vice-President

中国科学院海西研究院转化医学中心

Translational Medicine Research Center, Haixi Institutes
Chinese Academy of Sciences

Report of COVID-19 Antigen Rapid Test Kit (Colloidal Gold) about Omicron variant strain

Date: 2021-12-02

Product Name: COVID-19 Antigen Rapid Test Kit (Colloidal Gold)

Package: 1 Test/Kit

REF. No.: CG01Ag-01

Manufacturer: Xiamen AmonMed Biotechnology Co., Ltd.

Company Address: Unit 503, No. 120 Xinyuan Road, Haicang District, Xiamen, Fujian, China

Translational Medicine Research Center, Haixi Institutes, Chinese Academy of Sciences

Authorized Signature & Seal:



Tel: +86-592-3594011/01

E-mail: xiaozhang@fjirms.ac.cn

<http://www.xmirem.ac.cn/kydw/kytd/ktz9/ktz9/>

No. 258 Duishan Xiheng Road, Jimei District, Xiamen 361021, P.R. China

中国科学院海西研究院转化医学中心

Translational Medicine Research Center, Haixi Institutes

Chinese Academy of Sciences

Information of AmonMed COVID 19 Antigen Rapid Test Kit (Colloidal Gold) about Omicron variant strain

The Translational Medicine Research Center was conducted from Alpha variant strains base on routine bioinformatics analysis and laboratory verification have been carried out. Computer stimulation data and key mutant sites pseudovirus experimental results have comfrimed that AmonMed COVID 19 Antigen Rapid Test Kit (Colloidal Gold) can detect major Novel coronavirus variant strain, including the Delta variant strain which has been validated by a large number of clinical tests. According to bioinformatics comparison of the novel Coronavirus "Omicron" variant strain found in South Africa at present. The new variant strain B.1.1.529 has about more than 50 mutation sites, including 32 mutation sites on S protein and 4 mutation sites on N protein, of which 3 mutation sites have appeared in the last year. Therefore, AmonMed COVID 19 Antigen Rapid Test Kit (Colloidal Gold) can detect Omicron mutant strain.

Conclusion

AmonMed COVID 19 Antigen Rapid Test Kit (Colloidal Gold) can detect Omicron mutant strain.



Tel: +86-592-3594011/01

E-mail: xiaozhang@fjirsm.ac.cn

<http://www.xmirem.ac.cn/kydw/kytd/ktz9/ktz9/>

No. 258 Duishan Xiheng Road, Jimei District, Xiamen 361021, P.R. China

Certificate of Registration



This is to certify that the Quality Management System of

Xiamen AmonMed Biotechnology Co., Ltd.

Unified Social Credit Code: 913502050899181912

Operation Address: Unit 503 & 1203, 120 Xinyuan Road, Haicang District, Xiamen City, Fujian Province, China(Production); 5F and 6F, No.253, Duiying South Road, Jimei District, Xiamen City, Fujian Province, China(Production, Office)

Registered Address: Unit 503, 120 Xinyuan Road, Haicang District, Xiamen City, Fujian Province, China

applicable to

Production and sales of in vitro diagnostic reagents (within the scope of qualification, see attachment for details); production and sales of COVID-19 IgM/IgG test kit(Colloidal Gold); COVID-19 antigen rapid test kit(Colloidal Gold), COVID-19 neutralizing antibody test kit(Colloidal Gold), COVID-19/influenza A /influenza B virus antigen assay kit(Colloidal Gold)(export to EU)

has been assessed and registered by NQA against the provisions of

ISO 13485:2016

This registration is subject to the company maintaining a quality management system, to the above standard, which will be monitored by NQA.

Certified Clients shall accept regular surveillance assessments, the validity of certificates shall be maintained for the positive result of audit.

The information of this certificate can be checked on CNCA's website (www.cnca.gov.cn)

SNQA's website: www.snqa.com.cn

Managing Director

Certificate Number

46652

Date:

09 July 2019

Reissue Date:

01 July 2021

Valid Until:

09 July 2022



0015



The use of the UKAS Accreditation Mark indicates accreditation in respect of those activities covered by the accreditation certificate number 015 held by NQA.

NQA is a trading name of NQA Certification Limited, Registration No 09351758. Registered Office: Warwick House, Houghton Hall Park, Houghton Regis, Dunstable, LU5 5ZX, UK.

This certificate is the property of NQA and must be returned on request.



Xiamen AmonMed Biotechnology Co., Ltd.

Annex

1. Kit for the Determination of D-Dimer (Fluorescence Immunochromatography)
2. Kit for the Determination of Procalcitonin (Fluorescence Immunochromatography)
3. Kit for the Determination of Urine Microalbumin(U-ALB) (Fluorescence Immunochromatography)
4. Kit for the Determination of The whole course of C-Reactive Protein (hsCRP+CRP) (Fluorescence Immunochromatography)
5. Bacterial Vaginosis Test Kit (Sialidase)
6. Occult Blood (Hemoglobin/Transferrin) Test Kit (Colloidal Gold)
7. Kit for the Determination of Troponin I (Fluorescence Immunochromatography)
8. Kit for the Determination of Credine Kinase-MB (Fluorescence Immunochromatography)
9. Kit for the Determination of Myoglobin (Fluorescence Immunochromatography)
10. Kit for the Determination of Cystatin C (Fluorescence Immunochromatography)
11. Kit for the Determination of N-terminal pro-brain natriuretic (Fluorescence Immunochromatography)
12. Kit for the Determination of Hemoglobin Alc (Fluorescence Immunochromatography)
13. Kit for the Determination of Heart-type Fatty Acid Binding Protein(Fluorescence Immunochromatography)
14. Kit for the Determination of $\beta 2$ Microglobulin (Fluorescence Immunochromatography)
15. Kit for the Determination of Neutrophil gelatinase-associated (Fluorescence Immunochromatography)
16. Kit for the Determination of 25-OH Vitamin D(Fluorescence Immunochromatography)
17. Kit for the Determination of Troponin I /CKMB/Myoglobin(Fluorescence Immunochromatography)
18. COVID-19 IgM/IgG Test Kit (Rare Earth Nano Fluorescence Immunochromatography)

Managing Director



0015



Certificate Number **46652**

Date: 09 July 2019
Reissue Date: 01 July 2021
Valid Until: 09 July 2022





COVID-19 Antigen Rapid Test Kit (Colloidal Gold)

Instructions for Use with Saliva Swab Specimen

For self-testing/home use/private use

【INTENDED USE】

This test kit is used for in vitro qualitative detection of SARS-CoV-2 antigens in human saliva swab samples. It is intended for rapid detection of suspected COVID-19 cases within the first 7 days of symptom onset.

A positive test result indicates that the sample contains SARS-CoV-2 antigen. A negative test result does not rule out the possibility of infection.

This test kit is for self-testing by lay person in a non-laboratory setting (such as user's home or certain non-traditional sites such as airports, offices, schools, stadiums, etc.). The test results of this test kit are for preliminary screening and clinical reference only. It is recommended to conduct a comprehensive analysis of the condition based on the user's clinical manifestations and other laboratory tests.

【TEST PRINCIPLE】

This kit uses immunochromatography for detection. The specimen will move forward along the test card under capillary action. If the specimen contains a novel corona virus antigen, the antigen will bind to the colloidal gold-labeled new corona virus monoclonal antibody. The immune complex will be captured by corona virus monoclonal antibodies which are membrane fixed, from the fuchsia line, display will be corona virus antigen positive, if the line does not show color, the negative result will be displayed. The test card also contains a quality control line C, which shall appear fuchsia regardless of whether there is a detection line.

【MATERIALS PROVIDED】

Components	Specification		
	1 Test/Kit	5 Tests/Kit	25 Tests/Kit
	CG01Ag-01S-ST	CG01Ag-05S-ST	CG01Ag-25S-ST
Test card	1	5	25
Extraction solution	1	5	25
Saliva swab	1	5	25
Extraction tube	1	5	25
Instructions for use	1	1	1
Tube rack	1(packaging)	1	1

【PERFORMANCE CHARACTERISTICS】

Sensitivity: 96.55% (95% CI, 93.05% -98.32%)

Sensitivity: The true positive rate
Specificity: >99% (95% CI, 99.19% -100.00%)
Specificity: The true negative rate
Accuracy: 98.86% (95% CI, 97.87% -99.50%)
Accuracy: The true negative and positive rate
Limit of Detection: 5×10²TCID₅₀/mL

【Cross-reactivity】
 The sample with human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, adenovirus, human metapneumovirus, parainfluenza virus 1, parainfluenza virus 2, parainfluenza virus 3, parainfluenza virus 4, and etc. have no cross reaction.

【INTERFERENCES】
 Common interfering substances in the sample, such as blood, mucin and etc. have no effect on the test results.

- 【WARNINGS AND PRECAUTIONS】**
- Children under 18 years of age should be assisted by an adult.
 - Read the Instructions for Use (this leaflet) carefully before use.
 - Do not re-use. Do not drink any liquid in the test kit.
 - Do not use the test kit beyond the expiry date.
 - Do not use the test kit if any of the kit components are missing, broken, or unsealed.
 - Store the test kit at 2-30°C. Do not freeze.
 - Handle all specimens as potentially infectious.
 - The specimens should be tested immediately after collection.
 - Inadequate or inappropriate specimen collection, storage, and transport may yield inaccurate test results.
 - Correct specimen collection is a quite important step during the testing procedures. Make sure to collect enough specimens with the saliva swab.
 - The test should be used at room temperature (8-30 °C). If the test has been stored in a cool area (less than 8 °C), leave it at normal room temperature for 30 minutes before using.
 - Use the saliva swab provided in the test kit to ensure optimal performance of the test.
 - Apply the drops of test specimen only to the specimen well (S) on the test card.
 - Too many or too few drops of extraction solution may result in invalid or incorrect test result.
 - The specimen collection procedures may be uncomfortable. Do not insert the saliva swab too much deeper, please stop the test if you feel strong resistance or pain.
 - Keep the test kit and kit components out of the reach of children and pets before and after use.
 - Wear safety mask or other face covering when collecting saliva swab specimen from child or another individual.
 - Use of gloves is recommended when conducting testing.

【WASTE DISPOSAL AFTER TEST PROCEDURES】

- Place the used test card, extraction solution and saliva swab in a disposal bag and seal the disposal bag.
- Dispose all used devices and other components into normal household waste container in compliance with the applicable local regulations.
- Wash or sanitize your hands again.

【INTERPRETATION OF TEST RESULT】

Positive:
 If both the control line (C) and the test line (T) appear within 15-20 minutes, the result is positive.

Caution: No matter how faint the colored band is in the test line(T), the result should be considered as positive.

Negative:
 If there is only a control line (C) and test line (T) is colorless within 15-20 minutes, the test result is negative.

Invalid:
 If the control line (C) is not observed within 15-20 minutes, the test is invalid. And the test shall be conducted again with a new test card.

【FREQUENTLY ASKED QUESTIONS (FAQ)】

- When can/should I test myself?
 You can have a test on yourself whether you have symptoms or not. Please note that the test result is a snapshot that is valid for this point in time. Tests should therefore be repeated according to local regulations.
- What should I pay attention to in order to have the optimal test result?
 Always follow the instructions for use correctly. Perform the test immediately after collecting the sample. Apply two drops from the extraction tube into the specimen well of the test card. Too many or too few drops can lead to an incorrect or invalid test result.
- The test strip is very discolored. What may be the reasons?
 The reason for a clearly visible discoloration of the test strip is that too many drops has been dispensed from

【LIMITATIONS】

- The components of this test kit are to be used exclusively for the qualitative detection of SARS-CoV-2 antigen in saliva swab specimens. Other specimen types may lead to incorrect results and must not be used.
- The test kit is used for rapid detection of suspected COVID-19 cases within the first 7 days of symptom onset, so asymptomatic individuals may get a false-negative test result.
- Failure to follow the instructions for test procedures and interpretation of test results may adversely affect test performance and/or produce invalid results.
- A negative test result may occur if the specimen was collected or extracted improperly. A negative test result does not eliminate the possibility of SARS-CoV-2 infection and should be confirmed by a molecular assay.
- Improper storage, collection, or even freezing and thawing of the specimen can lead to inaccurate test results.
- Positive test results do not rule out co-infections with other pathogens.
- If the viral load of the specimen is below the detection limit of the test, the test may produce a negative result.
- Test results must be evaluated in conjunction with other clinical data available to the physician laboratory test results.
- The amount of antigen in a sample may decrease as the duration of illness develops. Specimens collected after 5-7 days of symptom onset of illness are more likely to be tested negative compared to a molecular assay.

【STORAGE AND SHELF LIFE】

- The test kit should be stored at 2-30°C, and the shelf life is 18 months.
- After the aluminum foil pouch is unsealed, it is recommended to use the test card within 1 hour at room temperature.
- The extraction solution is recommended to be used within 1 hour after opening at room temperature.

【PREPARATION BEFORE TEST PROCEDURES】

- Make sure all kit components are equilibrated to room temperature on the flat and clean surface.
- Make sure the kit components are complete without any missing or damaged after opening.
- Make sure to check the kit expiry date before testing.
- Make sure to wash or sanitize your hands, and make sure they are dry before starting.
- Make sure to prepare the following materials required but not provided in the kit.
 - Timer (watch)
 - Any necessary personal protective equipment (gloves, glasses etc.)
 - Waste container

【OPERATION OF TEST PROCEDURES】

- Take out the Instructions for Use and read it carefully.

- the extraction tube into the specimen well of test card. The indicator strip can only hold a limited amount of liquid. If the control line (C) does not appear or the test strip is very discolored, please retest by using a new test card according to the instructions for use.
- I have taken the test, but the control line (C) doesn't appear. What should I do?
 According to the instructions for use, this test result is invalid. Please retest by using a new test card.
 - I am not sure about reading test result. What should I do?
 Read the instructions for use again, and if this doesn't help, please contact the nearest health facility recommended by your local authorities for help.
 - If my test result is positive, what should I do?

- Take out the tube rack and assemble it. Gently press one tube rack well and place the extraction tube into the tube rack.
 Note: For specification of 1 Test/Kit, tube rack is on the kit packaging

- Hold the extraction solution vial with your fingers and make sure the tail is upward. Rotate the tail of the extraction solution vial.
 Caution: Safely unscrew the vial away from your eyes and face. Be careful of the sharp edge of the vial. Do not pour out the liquid.

- Squeeze all extraction solution from the vial into the extraction tube.
 Caution: Avoid touching the vial against the tube.

- Find the saliva swab in the sealed wrapper. Identify the fabric, soft tip of the saliva swab. Peel off the swab packaging and gently take out the saliva swab.
 Caution: Never touch the fabric, soft tip of the saliva swab with your fingers to avoid pollution.

- Specimen Collection
 Do not eat or drink anything, such as gum, tobacco, liquor, etc. 30 minutes prior to sampling.
- Insert the saliva swab by one hand into the mouth cavity.

- Place the saliva swab tip between upper and lower molar teeth, then gently bite the swab tip with upper and lower molar teeth for no less than 10 seconds and meanwhile close the mouth for complete saliva absorption in the depths of the mouth.
- After saliva collection, gently take out the swab.

NOTE: False negative results may occur if the saliva

There is possibility of hospitalization, complications and even death after infection with SARS-CoV-19. You should immediately contact the nearest health facility recommended by your local authorities.


- If my test result is negative, what should I do?
 If you test result is negative by the test, you also need to obey the local regulations. If you experience symptoms such as fever, headaches, migraines, loss of sense of smell and taste, contact the nearest health facility recommended by your local authorities.
- Will this test hurt?
 No, the saliva swab is not sharp and it should not hurt. Sometimes the saliva swab can make slightly uncomfortable or tickly. If you feel pain, please stop the test and ask for help from a healthcare provider.

specimen is not collected properly.



7. Specimen Handling

- Insert the saliva swab into extraction tube. Stir the saliva swab more than 5 times. Leave saliva swab in extraction tube for about 1 minute.
- Squeeze the swab against the inner wall of extraction tube to release the liquid as much as possible when you remove the swab. Dispose of the test swab with normal household waste in accordance with applicable local regulations.



8. Press the cap onto the extraction tube tightly.



9. Unseal the foil pouch and take out the test card. Place the card on the flat surface.



10. Apply 2 drops of extracted specimens to the specimen well of the test card, and then start timing.



11. Read the test results in 15-20 minutes, and test results after 20 minutes may not be accurate.



【ACCESSORY】

Accessory	Manufacturer	EC-Representative	CE-Mark
Saliva Swab	Shenzhen Kangdaan Biological Technology Co. Ltd. 3rd floor, Building A2, Shunheda factory, Liuxiandong industrial zone, Xilli street, Nanshan district, Shenzhen, China.	Share Info Consultant Service LLC Repräsentanzbüro Heerder Lohweg 83 40549 Düsseldorf, Deutschland	acc. 93/42/EEC

【EXPLANATION FOR SYMBOLS】

	Expiry date		Batch Number		See Instructions for use
	Test (s) per kit		Store at 2-30°C		Catalogue Number
	Manufacturer		CE Mark		Do not reuse
	In Vitro diagnostic use		European Authorized Representative		

【ISSUE DATE AND VERSION NO.】

Issue Date: Oct 15th, 2021; Version 4.0

Xiamen AmonMed Biotechnology Co., Ltd.
 Address: Unit 503, 120 Xinyuan Road, Haicang District, Xiamen, Fujian, China.

EC REP **SUNGO Europe B.V.**
 Address: Olympisch Stadion 24,1076DE Amsterdam, Netherlands