No. IFU-COVIDAg-NN-01, Ver. A/0

COVID-19 (SARS-CoV-2) Antigentestkit (kolloidales Gold)

REF SL030101

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[Microbial Interference Studies]

Microbial interference in the DeepBlue SARS-CoV-2 Ag Test was evaluated by testing a panel of related pathogens, high prevalence disease agents and normal or pathogenic flora to demonstrate that false negatives do not occur when SARS-CoV-2 is present in a specimen with other microorganisms.

Mikroorganismus	Konzentration	Interferenz (Ja / Nein)
Parainfluenzavirus Typ 1	1 x 10 ⁵ PFU/mL	Nein (3/3 positiv)
Parainfluenzavirus Typ 2	1 x 10 ⁵ PFU/mL	Nein (3/3 positiv)
Parainfluenzavirus Typ 3	1 x 10 ⁵ PFU/mL	Nein (3/3 positiv)
Parainfluenza-Virus Typ 4a	1 x 10 ⁵ PFU/mL	Nein (3/3 positiv)
Adenovirus (z. B. C1 Ad. 71)	1 x 10 ⁵ PFU/mL	Nein (3/3 positiv)
Humanes Metapneumovirus (hMPV)	1 x 10 ⁵ PFU/mL	Nein (3/3 positiv)
Influenza A H3N2(Wisconsin/67/05)	8.82 x 10 ⁴ PFU/mL	Nein (3/3 positiv)
Influenza A H1N1	1 x 10 ⁵ PFU/mL	Nein (3/3 positiv)
Haemophilus influenzae	1 x 10 ⁶ CFU/mL	No(3/3 positive)

Streptococcus pneumoniae	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Streptococcus pyogenes	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Influenza B (Malaysia/2506/04)	2.92 x 10 ⁴ PFU/mL	No (19/20 positive)
Enterovirus	1 x 10 ⁵ PFU/mL	No (3/3 positive)
Respiratory syncytial virus	1 x 10 ⁵ PFU/mL	No (3/3 positive)
Rhinovirus	4.17 x 10 ⁵ PFU/mL	No (3/3 positive)
Chlamydia pneumoniae	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Legionella pneumophila	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Mycobacterium tuberculosis	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Pneumocystis jirovecii	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Pseudomonas Aeruginosa	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Candida albicans	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Pooled human nasal wash	14% v/v	No (3/3 positive)
Bordetella pertussis	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Mycoplasma pneumoniae	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Staphylococcus Epidermidis	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Streptococcus Salivarius	1 x 10 ⁶ CFU/mL	No (3/3 positive)
Human coronavirus 229E	1 x 10 ⁵ PFU/mL	No(3/3 positive)
Human coronavirus OC43	1 x 10 ⁵ PFU/mL	No (19/20 positive)
Human coronavirus NL63	9.87 x 10 ³ PFU/mL	No(3/3 positive)
MERS coronavirus	7930 PFU/mL	No (3/3 positive)

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Reference RT-PCR Assay						95% Wilson Score CI			
						LCI	UCI		
DEEP BLUE SARS-CoV-2 Ag Test		POS	NEG	Total	PPA	96.4%	90.8%	98.2%	
		POS	106	1	107	NPA	99.8%	94.4%	99.9%
		NEG	4	409	413	PPV	99.1%	93.7%	99.8%
		TOTAL	110	410	520	NPV	99.0%	93.5%	99.7%

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LCI – unteres Konfidenzintervall

UCI – oberes Konfidenzintervall

ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO.,LTD.
 4th Floor,D-1# Zone, Pearl Industrial Park, 106 Innovation Avenue,High-Tech Development Zone,230088 Hefei, Anhui,China

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Anhui Deepblue Medical Technology Co., Ltd.

Website: www.dbluemedical.com

Address: 4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue, High-Tech Development Zone 230088 Hefei, Anhui, China

COVID-19 (SARS-CoV-2) Antigen Test kit (Colloidal Gold)

Verification reports of the clinical samples

Anhui DeepBlue Medical Technology Co., Ltd.

Web: www.dbluemedical.com

**Address: Innovation Avenue 106, High-Tech Development Zone
of Anhui Province 4th Floor D-1# Zone, Pearl Industrial Park,
230088 Hefei, Anhui, China**

Anhui Deepblue Medical Technology Co., Ltd.

Website: www.dbluemedical.com

Address: 4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue, High-Tech Development Zone 230088 Hefei, Anhui, China

1. Objections

- The clinical verification research of the COVID-19 (SARS-CoV-2) Antigen Test kit is performed by comparing the Antigen test results with the commercialized and authorized comparison kit test results, so as to confirm the clinical application performance of the antigen test kit. Based on the peculiarity of the SARS-CoV-2 and the related requirements of CE and WHO, the collected samples were prospectively tested and two kinds of test methods has been adopted for the determination of the same sample. The test results of both methods will be compared to verify the sensitivity and specificity of the Antigen test kit. And all test results are statistically analyzed.
- In the verification research, all the adopted positive samples have been tested and confirmed as positive by the real time fluorescent PCR kit. In the whole verification research, 110 positive samples and 410 negative samples have been adopted for the verification.

2. Methods

- Specimen category: nasopharyngeal swab
- 110 positive samples and 410 negative samples were tested.
- All nasopharyngeal swab samples were collected from the subjects who have been diagnosed with real time fluorescent PCR kits.
- All nasopharyngeal swab samples were tested and all results were recorded. After the test, all results were compared with those of the fluorescent PCR.
- **Design of the experiments process:** All the collected positive and negative samples will be randomly shuffled together and re-numbered. Then these re-numbered samples are sent to the test operators. Of note, it is ensured that the test operators only know the number of the samples and nothing about the detailed information about the samples. The test operators will test all adopted samples with our antigen test kit and record the No. and test results of each sample on the sample transfer sheet. After the completion of all tests, all sample transfer sheets will be collected. Based on the random number and the corresponding sample number of each

Anhui Deepblue Medical Technology Co., Ltd.

Website: www.dbluemedical.com

Address: 4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue, High-Tech Development Zone 230088 Hefei, Anhui, China

specimens, the test results of all double-blind samples will be unblinding. All verification results will be summarized and the test results of our antigen test kit will be compared to the control test kit to confirm the consistency between these two test kits.

- Test site and time: The Center of Disease Control and Prevention of Anhui Province (Anhui CDC)
- Time: From 09/10/2020 to 12/09/2020.

3. Sample collection

- The storage, transfer and using of storage buffer of the nasopharyngeal swab samples were all conducted according to the IFU of the test kit.
- The swab samples will be treated before the test according to the detailed operation procedures in the IFU of the test kit.

4. Principle of the test

- The principle of this antigen detection is based on the dual antibody sandwich test mechanism. Typically, the antibody against the N protein of the SARS-CoV-2 is labeled with the colloidal gold particles while the coating antibody against the N protein is sprayed onto the nitrocellulose membrane of the strip as the test line. These two antibodies can both recognize and bind with the N protein. The tested virus is of N protein and S protein. The structure of the N protein is more integrate than S protein. Detection of N protein can effectively judge the existence of the SARS-CoV-2, providing the direct proof for the diagnosis of the subject infected by the SARS-CoV-2 or not.
- All the samples have been confirmed with the SARS-CoV-2 nucleic acid detection kit (real time fluorescent PCR method) manufactured by the Sun Yet-Sen University DAAN Gene Co. Ltd.
- If the abnormal results occurred, the samples will be confirmed with the third-part real time fluorescent PCR kit of BGI (Wuhan) (authorized by the FDA of EUA of US and NMPA of China).

Anhui Deepblue Medical Technology Co., Ltd.

Website: www.dbluemedical.com

Address: 4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue, High-Tech Development Zone 230088 Hefei, Anhui, China

- Test kit: COVID-19 (SARS-CoV-2) Antigen Test Kit (Colloidal Gold)
- Manufacturer: Anhui DeepBlue Medical Technology Co. Ltd.
- Size: 25 pc/box
- Production Batch No.: 200901

5. Incompatibility analysis

- If there were some inconsistent results obtained in the clinical research, the inconsistent test plan should be made and conducted.
- For the results comparison study with the real time PCR method, the results of two methods should be consistent. If the results were inconsistent, the third-part real time fluorescent PCR kit for BGI (Wuhan) Co. Ltd. (authorized by FDA of US and NMPA of China) should be adopted for the confirmation.
- Variation analysis should not be used for the change of the performance data and it can be adopted and added as the footnote of the performance table.
- If necessary, for the confirmation of the sensitivity, the positive samples can be diluted serially to demonstrate the analysis sensitivity is decreased or increased with the comparison protocol.

6. Statistics of the results

- Statistics of the patient information (gender, age, symptoms, Ct of PCR) for the verification research of 520 samples.

The table shows the positive samples categorized according to the onset time of the patients:

Onset time (day)	DeepBlue COVID-19 Ag (n=110)		
	Amount of Sample	Antigen positive	Sample percentage
0-3 days	44	40	40%
4-7 days	44	43	40%
> 7 days	22	20	20%

Anhui Deepblue Medical Technology Co., Ltd.

Website: www.dbluemedical.com

Address: 4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue, High-Tech Development Zone 230088 Hefei, Anhui, China

A. Clinical performance evaluation standards

- Compared with the commercial nucleic acid detection kits (verified), to verify the positive percentage agreement (PPA) and negative percentage agreement (NPA) of our antigen test kit in the clinical applications. The detection rate of our antigen test kit is compared with commercial nucleic acid detection kit.
- The verification kit and the control kit were all conducted based on the recommended operation procedure in the operation manual in the kit.

B. Statics of the clinical test results

- Sensitivity (positive detection percentage)

Table1. Test results of 110 positive nasopharyngeal samples

Test results	Antigen detection	RT-PCR
Positive (+)	106	110
Negative (-)	4	0
Total	110	110
Accuracy (Sensitivity)	96.4%	100%

- Specificity (negative detection percentage)

Table 2. The test results of 410 negative nasopharyngeal samples

Test results	Antigen detection	RT-PCR
Positive (+)	1	0
Negative (-)	409	410
Total	410	410
Accuracy (Specificity)	99.8%	100%

- In the above verification study, 100 positive and 410 negative swab samples were collected, and the whole sample volume is 520 cases. Among the positive samples, the Ct values of 20 samples were high than 30 (>30).
- For the 5 samples with abnormal results, these 5 samples were re-tested with the third-part real time PCR kit and the results were in consistent with the original test results. The reason for this is due to the difference

Anhui Deepblue Medical Technology Co., Ltd.

Website: www.dbluemedical.com

Address: 4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue, High-Tech Development Zone 230088 Hefei, Anhui, China

between the DeepBlue antigen test kit and the reference detection kit.

With the statistic analysis, there was no significant variations.

Table 3. Statistical analysis of the test results of the nasopharyngeal swab samples

Reference RT-PCR							95% CI	
							LCI	UCI
DeepBlue COVID- 19 Ag test		POS	NEG	Total	PPA	96.4%	90.8%	98.2%
	POS	106	1	107	NPA	99.8%	94.4%	99.9%
	NEG	4	409	413	PPV	99.1%	93.7%	99.8%
	Total	110	410	520	NPV	99.0%	93.5%	99.7%

PPA-positive percentage agreement (sensitivity)

NPA-negative percentage agreement (specificity)

PPV-positive predictive value

NPV-negative predictive value

CI- Confidential interval

LCI-low confidential interval

UCI-up confidential interval

C. Test results

Test results demonstrate that: In this verification research, 110 positive and 410 negative nasopharyngeal swab samples were tested respectively. The **sensitivity is 96.4%** and the **specificity is 99.8%**. The test results of nasopharyngeal swab samples are basically in consistent with those of the nasopharyngeal swab samples.