

Date: March 31, 2023

## Statement

To whom it may concern,

**We, Beijing Wantai Biological Pharmacy Enterprise Co., Ltd.** registered at **No. 31 Kexueyuan Road, Changping District, Beijing 102206, China**, as one of the leading Chinese enterprises specialized great progress in both product development and market presence. We've been committed to providing professional IVD solution to customers over the world.

All viruses, including SARS-CoV-2, the virus that cause COVID-19, change over time. Some changes may effect the virus's properties, such as how easily it spreads, the associated disease severity, or the performance of vaccines, therapeutic medicines, diagnostic tools, or other public health and social measures. The identified SARS-CoV-2 variants are as follows:

WHO label	Pango lineage	GISAID clade/lineage	Nextstrain clade
Alpha	B.1.1.7	GRY	20I(V1)
Beta	B.1.351	GH/501Y.V2	20H(V2)
Gamma	P.1	GR/501Y.V3	20J(V3)
Delta	B.1.617.2	G/478K.V1	21A
Lambda	C.37	GR/452Q.V1	21G
Eta	B.1.525	G/484K.V3	21D
Iota	B.1.526	GH/253G.V1	21F
Kappa	B.1.617.1	G/452R.V3	21B
Mu	B.1.621	GH	21H
Omicron*	B.1.1.529	GR/484A	21K, 21L, 21M, 22A, 22B, 22C

\*Includes BA.1, BA.2, BA.3, BA.4, BA.5 and descendent lineages such as BA.2.12, BA.2.12.1, BA.2.75, BA.5.2, BA 2.75.2, BQ.1, BQ 1.1, XBB.1, XBB.1.5, XBB.1.16, etc.

hereby we declare that, the above-listed variants can be tested by WANTAI SARS-CoV-2 Ag Rapid Test (Colloidal Gold), code WJ-2901, WJ-2905, WJ-2910, WJ-2925 and WJ-2950, and WANTAI SARS-CoV-2 RT-PCR Rapid test, code WS-2248.

 Sincerely,  
Beijing Wantai Biological Pharmacy Enterprise Co., Ltd.