



SHENZHEN MICROPROFIT BIOTECH CO.,LTD
Room 1001 of Unit 2, Room 1001 and Room 1101 of Unit 1,
Building 2, Hongchuang Technology Center, Xikeng Community,
Fucheng Sub-district, Longhua District, 518000 Shenzhen,
Guangdong, PEOPLE'S REPUBLIC OF CHINA
Tel: (+86)755 29303416 Fax: (+86)755 61688111
Http://www.microprofit-bio.com

Statement

To Whom it might concern

We, Shenzhen Microprofit Biotech Co., Ltd., as the manufacturer of SARS-CoV-2 mutation Impact on SARS-CoV-2 & Influenza A/B Antigen Combo Test Kit(Colloidal Gold Chromatographic Immunoassay) (fluorecare - REF MF-69), SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit (Colloidal Gold Chromatographic Immunoassay) (fluorecare - REF MF-71),SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit (Colloidal Gold Chromatographic Immunoassay) (Self-test) (fluorecare - REF MF-71) hereby declare.

1.Further In-vitro studies to investigate the impact of mutated SARS-CoV-2 N proteins and the analysis performance of the product is finished, B.1.1.7, B.1.351, B.1.429, B.1.427, B.1.2, P.1, P.2, P.3, B.1.617.2,B.1.617.3, C.37, B.1.617.1, B.1.525, B.1.526, B.1.526.1, B.1.526.2, B.1.1.529, BA.1, BA.2, XE, BA.4, BA.5,BA.5.2, BA.4.6, BF.7, BA.2.75, B.1.1.318, C.36.3, AY.4.2, BA.2.10, XBB, XBB.1.5, XBC, BQ.1, BQ.1.1, BJ.1,BA.2.3.20, B.1.621, BM.2, BA.4.7, BN.1, BA.5.2.35, BQ.1.4, BQ.1.3, BQ.1.2, BA.5.9, BM.1.1.1 BS.1,BA.2.75.5,BA.2.75.2, BA.2.12.1, BA.2.86, EG.5.1, JN.1,KP.2,KP.3.1.1,XEC,NB.1.8.1 the N protein mutation site is not in the recognition region of the coated and labeled antibody, can recognize the reaction well.

2. In the future, our company will continue to follow up the mutation of the new coronavirus, and timely evaluate and verify the detection ability of mutated recombinant protein and clinical performance to ensure the sensitivity and specificity of the detection kit.

Shenzhen Microprofit Biotech Co., Ltd

General Manager: Ying Liu

Signature:

Date:

7075.0528

Ying Liu